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21011 PASCAL ball

1. Purpose

The product is mainly for the high school physics demonstration liquid transfer to the pressure.

2. Structure.

At one end of the aluminum alloy cylinder, mounded with a hollow metal ball, ten steel gas nozzle on hollow ball, with the piston within the cylinder, piston by cowhide bowl bowl (or oil resistant rubber rings). And move up and down inside the cylinder, piston on rod out cover and equipped with a handle.

3. Usage

Unscrew the hollow ball, the cylinder into the water, bring the piston, the cylinder under the suction, the hollow ball is filled with water, and put the ball in the water cycle to the cylinder, then the proposed the water slowly down the piston, water injection in all directions from the nozzle holes, showed that the liquid have received the pressure transmission in all directions.

Note: while doing the experiment to the line of sight arranged nozzle jets direction, avoid the water spray teachers themselves and students.

4. Maintain

Finished the experiment instruments should be apart rub-up, air drying after assembly. Metal parts with a few vaseline, prevent corrosion.